

Interface and command description for the Control and Supply Unit 2 (CSU2) of IFG Institute for Scientific Instruments GmbH (as of 2016-03-21, translated July 2016)

General remarks regarding the network interface

Communication with the CSU2 works over a TCP/IP network connection. The CSU2 acts as TCP server using a fixed IP address and port (standard is 192.168.1.3 port 23, i.e. telnet port). In case a different IP address is needed, preferably inform IFG before delivery of the hardware.

Overview of command syntax

Commands start with a dollar sign '\$', immediately followed by a command mnemonic (two upper case letters, digits, or special (printable) characters). One or more parameters may follow. Each command ends with a carriage return character (in the following represented as <CR>). All characters between two commands ('<CR>...\$') are ignored, i.e. commands can also be completed with <CR><LF> (carriage return, line feed).

Parameter can be whole numbers (no float values), switches (boolean values) or other characters. Switches can be 'ON' and 'NO' or alternatively '+' and '-'. Both forms are equivalent. The first parameter may act as a subcommand, for example see the HV group of commands.

A parameter may be preceded with blanks. Two parameters must be separated with at least one blank. The last parameter is directly followed by <CR>.

Every '\$' character starts a new command, i.e. previously read characters get discarded silently. A command must not contain more than 89 characters.

There is a **response** to each command. The response might be an error message or a command confirmation. A response starts with an exclamation mark '!' instead of a '\$' sign.

A list of commands follows. Parameters are represented as lower case strings. Optional parameters are bracketed in '[...]'.

The query command **\$OK** returns the time since power-on (in ms). The value returned is a double word (4 bytes) with a maximum of 4 294 967 295. After that it starts again from zero. Its primarily purpose it to check synchronicity of computer and CSU2.

For **remote control** using the interface described here the CSU2 has to be in "PC" mode: The LED next to the green button "PC" should be illuminated. The actual mode can be queried using the command **\$RM?**. **All modifying commands lead to a positive response in local mode „LOC“ – yet they are not executed!** Read only commands function in local mode just as in remote mode.

To ensure an adequate **warm-up** of the tube and operation in line with its **prefomance limits** the type of tube (see TT-group of commands), and system time (\$RKOK, \$RKTP) have to be set correctly.

Command syntax

1	OK	14	FH	27	RK TP
2	HV	15	FU	28	RK OK
3	HV W-	16	XR	29	RK R?
4	HV WP	17	XR ?	30	RK T?
5	HV W?	18	RM ?	31	RK L?
6	HV UP	19	TT I?	32	RK S?
7	HV IP	20	TT IP	33	RK RX
8	HV HP	21	TT #?	34	RK RN
9	HV U?	22	TT LP	35	RK CX
10	HV I?	23	RK PP	36	RK CN
11	HV H?	24	RK LP	37	RK ID
12	HV ??	25	RK PB	38	XV
13	HV ?1	26	RK LB	39	<<

1 Command OK: Query system time since power-on

Query: \$OK<CR>
Response: !OK system time in ms<CR>

2 Command HV: Power high voltage generator on/off

Query: \$HV cmd<CR> cmd = ON bzw. '+': on
Response: !HV<CR> NO bzw. '-': off

3 Command HVW-: Skip warm-up of tube (immediately "warmed-up")

Query: \$HVW-<CR>
Response: !HVW-<CR>

4 Command HVWP: Manually set last time the tube was used

Query: \$HVWP cmd<CR> cmd = Time in format YYYY-MM-DD-hh:mm:ss
Response: !HVWP<CR>

- IfG - Institute for Scientific Instruments GmbH • Rudower Chaussee 29/31 • 12489 Berlin • Germany

17 Command XR?: Query shutter state

```
-----
Query: $XR?<CR>                status = ON bzw. '+': shutter is open
Response: !XR status<CR>        NO bzw. '-': shutter is closed
```

18 Command RM?: Query remote mode

```
-----
Query: $RM?<CR>                status = ON bzw. '+': remote control on
Response: !RM status<CR>        NO bzw. '-': local mode
```

19 Command TTI?: Query index of selected tube type

```
-----
Query: $TTI?<CR>                ix = Index 0..(ixcnt-1)
Response: !TTI? ix<CR>
```

20 Command TTIP: Set tube type by index

```
-----
Query: $TTIP ix<CR>            ix = Index 0..(ixcnt-1)
Response: !TTIP<CR>
```

21 Command TT#?: Query number of known tube types

```
-----
Query: $TT#?<CR>                ixcnt = biggest index + 1
Response: !TT#? ixcnt<CR>
```

22 Command TTLP: Query caption of tube type at index ix

```
-----
Query: $TTLP ix<CR>            ix = index 0..(ixcnt-1)
Response: !TTLP cap<CR>        cap = caption/description
```

23 Command RKPP: Program a string

```
-----
Query: $RKPP position string<CR>    position = 0..99 (101)
Response: !RKPP<CR>                String length of max. 32 Byte
```

24 Command RKLP: Read string from position

```
-----
Query: $RKLP position <CR>          position = 0..101
Response: !RKLP position String<CR> string length of max. 32 Byte
```

25 Command RKPB: Program pair of words

```
-----
Query: $RKPB position value1 value2<CR>    position = 0..99
Response: !RKPB<CR>                value = 0..65535
```

26 Command RKL B: Read pair of words from position

```
-----
Query: $RKL B position<CR>          position = 0..99
Response: !RKL B position value1 value2<CR> value = word (2 bytes)
```

27 Command RKTP: Set realtime clock

```
-----
Query: $RKTP cmd<CR>              cmd = time as YYYY-MM-DD-hh:mm:ss
Response: !RKTP<CR>
```

28 Command RKOK: Query realtime clock

```
-----
Query: $RKOK<CR>
Response: !RKOK cmd <CR>          cmd = time as YYYY-MM-DD-hh:mm:ss
```

```

29 Command RKR?: Query tube temperature
-----
Query: $RKR?<CR>
Response: !RKR? temp<CR>                                temp = temperature in m°C

30 Command RKT?: Query temperature of HV generator (SIC)
-----
Query: $RKT?<CR>
Response: !RKT? temp<CR>                                temp = temperature in m°C

31 Command RKL?: Query temperature of LED board (PT100)
-----
Query: $RKL?<CR>
Response: !RKL? temp<CR>                                temp = temperature in m°C

32 Command RKS?: Query temperature of shutter board (PT100)
-----
Query: $RKS?<CR>
Response: !RKS? temp<CR>                                temp = temperature in m°C

33 Command RKRX: Query maximum tube temperature
-----
Query: $RKRX<CR>
Response: !RKRX temp<CR>                                temp = temperature in m°C

34 Command RKRN: Query minimum tube temperature
-----
Query: $RKRN<CR>
Response: !RKRN temp<CR>                                temp = temperature in m°C

35 Command RKCX: Query maximum temperature of HV generator
-----
Query: $RKCX<CR>
Response: !RKCX temp<CR>                                temp = temperature in m°C

36 Command RKC�: Query minimum temperature of HV generator
-----
Query: $RKC�<CR>
Response: !RKC� temp<CR>                                temp = temperature in m°C

37 Command RKID: Query serial numbers / seconds of operation
-----
Query: $RKID<CR>
Response: !RKID HVSerNo HVSeconds TubeSerNo TubeSeconds<CR>
                                   SerNo:  serial number (string[32])
                                   Seconds: seconds of operation

38 Command XV: query firmware version
-----
Aufruf : $XV <CR>                                sig : signature (string[16])
Response: !XVsig<CR>

39 Command <<: reset programmed parameters
-----
Query: $<<<CR>
Response: !<<

40 Command R#0: Query parameter version (number and type)
-----
Query: $R#0<CR>
Response: !R#9999

```

Error messages

!ERROR: 00 General error or buffer overflow
!ERROR: 01 Unknown command
!ERROR: 02 Numeric parameter expected
!ERROR: 04 Boolean parameter expected
!ERROR: 05 Additional parameter expected
!ERROR: 06 Unexpected parameter or character
!ERROR: 07 Illegal numeric value
!ERROR: 08 Unknown subcommand
!ERROR: 09 Function not implemented or no hardware support
!ERROR: 10 Flash-EPROM programming fault
!ERROR: 11 Error clearing Flash-EPROM
!ERROR: 12 Flash-EPROM read error
!ERROR: 13 Hardware error

Error codes (see command No. 13 \$HV?1)

Error codes are shown on the display “µA”. LEDs of “PC” and “LOC” buttons are constantly flashing.

Warnings:

1111 : Safety line at housing not connected (RJ45 connector)
1112 : External interlock (RJ45 connector at tube)
1113 : Interlock of HV-Generator
2111 : Real time clock broken or battery empty
2112 : Temperature of LED-board critical

Errors with self inhibition. These require power cycling the device:

3111 : Temperature sensor at LED-board broken / not connected
3112 : Temperature sensor at shutter-board broken / not connected
3121 : Temperature of LED-board above limit
3122 : Temperature of shutter-board above limit

3211 : HV LED at tube housing broken
3221 : Shutter hangs or shutter light bulb broken
3222 : Shutter hangs or shutter-LEDs at tube housing broken

3321 : Vacuum switch 1 broken
3322 : Vacuum switch 2 broken
3331 : HV powered on, and filament cable not (properly) connected
3332 : "PC"-mode active and HV on / shutter opened, but communication with PC timed out
3333 : No connection to HV generator. Power failure? Safety relais?